MOULDER[®]

KISS[®]

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Instructions for use

MOLDERM IS A MONOPHASIC, STERILE, NON-PYROGENIC AND INJECTABLE HYDROGEL FOR THE CORRECTION OF SOFT TISSUE DEFECTS. THE PRODUCT DEGRADES SLOWLY AFTER IMPLANTATION. THE HYDROGEL IS BASED ON HIGHLY PURIFIED CROSS-LINKED SODIUM HYALURONATE OF NON-ANIMAL ORIGIN, PRODUCED BY BACTERIAL FERMENTATION.

COMPOSITION

MOLDERM is available in three different concentrations of hyaluronic acid: 22 mg/ml, 24 mg/ml and 26 mg/ml.

MOLDERM TT in 1ml syringe, contains 22 mg/ml of cross-linked sodium hyaluronate, in sterile buffered water, with the addition of glycine and proline amino acids.

MOLDERM LIFT in 1 ml syringe, contains 24 mg/ml of cross-linked sodium hyaluronate, in sterile buffered water, with the addition of glycine and proline amino acids.

MOLDERM VOLUME in 1 ml syringe, contains 26 mg/ml of cross-linked sodium hyaluronate, in sterile buffered water, with the addition of glycine and proline amino acids.

MOLDERM KISS® in 1 ml syringe, contains 24 mg/ml of cross-linked sodium hyaluronate, in sterile buffered water, with the addition of glycine and proline amino acids.

Package contents

- 2 sealed blister packs, each containing 1 sterile single-use pre-filled syringe with luer-lock connection and 27G-13 mm or 25G-16 mm sterile needles.
- Patient information leaflet.
- Adhesive labels for traceability.
- Implant cards - Instructions for implant card use and correct filling in.

Indications

MOLDERM in the three different concentrations of hyaluronic acid is indicated for the restoration of facial volumes, correction of facial asymmetry, lipotrophy and scars.

Contraindications

MOLDERM must not be used in the following cases:

- In association with chemical peels, laser treatments or dermabrasion. The doctor will determine the waiting time necessary before giving the injection.
- In blood vessels
- In patients with a tendency to develop hypertrophic scars. In patients with infections and inflammations (acne, herpes, dermatitis, etc.) in situ or near the area to be treated.
- Patients treated with anticoagulants.
- Patients undergoing radiotherapy or ultrasound therapy in the area to be treated.
- Patients with hypersensitivity to hyaluronic acid.
- Patients with untreated epilepsy.
- Patients with known hypersensitivity or allergy to any of the components.
- Patients with a skin disease or alteration of any kind.
- Pregnant or breastfeeding women or on children

Precautions for use

The medical device is intended to be administered by qualified medical personnel by means of intradermal injection.

It is strongly recommended to obtain the patient's "informed consent" in advance.

Carefully check the patient's medical history before starting treatment. Explain the indications and the results that can be achieved to the patient.

Explain the contraindications and potential side effects of the treatment and advise on what to do in the immediate post-treatment period.

Respect the principles of asepsis and handling related to this type of procedure.

Write down the batch number of the product used on the traceability label, which will be attached to the patient's medical record.

It is advisable for the patient to avoid the intake of substances that affect the fluidity of the blood (aspirin, NSAIDs, Vitamin C and Vitamin E) near to treatment, in order to minimise the possibility of bruising or bleeding of the injected areas.

The use of the product should be absolutely excluded in areas where breast, tendon, bone or muscle implants are present.

For patients with a history or a declared autoimmune disease, the doctor will have to decide on the indication on a case-by-case basis, depending on the nature of the disease and the associated treatment. The doctor will also have to ensure careful monitoring of these patients, in particular by offering them a double test beforehand, and not injecting them if the disease is progressive.

Patients with a history of streptococcal disease, such as recurrent angina or acute joint rheumatism, should undergo a preliminary double test before any injection. In the case of acute joint rheumatism with cardiac localisation, it is recommended that the patient should not be injected.

Warnings

Do not use the product after its expiration date.

Do not use the product if the packaging is damaged, opened or stored improperly.

The contents of the pack are for single-use only and cannot be re-sterilised.

Once the package has been opened, the product must be used immediately: any remaining product must be disposed of.

Do not inject the product contained in the same syringe into more than one patient so as to avoid any risk of cross-contamination.

Do not mix the product with other substances.

Do not allow the product to come into contact with surgical instruments disinfected with quaternary ammonium bases such as benzalkonium chloride.

Do not add drugs to the product.

The contents of the syringe are disposable and its reuse can result in the loss of sterility with serious consequences for the patient's health. After use, syringes and needles can pose a potential biological hazard. Handle with care and dispose of in accordance with approved medical practice.

Dosage

Dosage depends on the condition of the patient's skin and the depth of the area to be treated. 0.5-1.0 ml per site are usually required. For each site, do not exceed the minimum quantity of a syringe (1 ml) per treatment session.

Directions for use

The treatment of adverse effects is the responsibility of the therapist. Before injecting the product, thoroughly disinfect the area to be treated with povidone-iodine or disinfectant containing chlorhexidine.

Unscrew the protective cap of the syringe.

Tighten the needle firmly onto the luer-lock connection of the syringe.

Gently press the plunger to remove any air from the product.

If the needle is obstructed, do not increase the pressure on the plunger, but replace the needle.

Slowly inject the product following the retrograde glide technique.

During the injection, the product can be shaped by the doctor with the contralateral hand.

When finished, massage the treated area gently to better distribute the product.

Post-treatment advice

The patient may be given an anti-inflammatory drug to reduce oedema (a physiological response to any implanted filler).

To reduce the inflammatory effects typically associated with an injection, ice can be applied to the treated region, even if not in direct contact with the skin.

Bruising/haematoma can be treated with appropriate creams or gels.

Do not press the skin on the treated area, nor apply cosmetics for 24 hours after applying the product.

Avoid surgery or dental procedures for two weeks after the date of injection. If the patient has to undergo these procedures, it is ALWAYS essential to inform the doctor concerned about the location of the implant.

Similarly, avoid beauty treatments (laser, peeling, etc.) in the treated area for 3 weeks from the date of implantation.

Avoid direct exposure to the sun and UVA-B rays (tanning lamps) for two weeks from the date of treatment.

Reactions to treatment and side effects

The doctor must inform the patient about the potential side effects resulting from the use of the product, which may occur immediately or after a certain period of time.

A number of adverse reactions after hyaluronic acid injection have been reported in the scientific literature. They include, but are not limited to: Inflammatory reactions (redness, oedema, erythema, etc.) possibly accompanied by itching or pressure pain or both;

- Bruising/haematoma;
- Indurations or nodules in the injection area;
- Stains and discolourations around the injection area;
- Intolerance reactions;
- Hypersensitivity in the treated area;
- Infection;
- Weak filling effect;
- Granuloma;
- Displacement of the material;
- Numbness;
- Lumps;
- Necrosis;
- Papules and nodules in the injection area;
- Allergic reactions

Immediately inform your doctor about any adverse event (mentioned above or otherwise) if it persists for more than 7 days. The implant must be checked and treated appropriately.

The resolution of any adverse event related to the implant is closely associated with the length of time from which the intervention occurred. Any side effects associated with treatment with MOLDERM must be reported to the manufacturer and/or distributor of the product.

Storage

Store the product at a temperature between 2°C and 30°C.

Do not expose the product to sunlight.

Do not freeze, do not heat.

The manufacturer accepts no responsibility for direct, indirect or consequential damages caused by the incorrect use or abnormal use of the medical device or by use in any manner other than that described in this patient information leaflet. The user is responsible for checking that the product is suitable for its intended use, and therefore assumes all risks and responsibilities in this regard.

Note: The gradation on the syringe has no measuring function but is given only as a visual check of the injected material.

Molderm Aesthetics AB
Dockplatsen 1, plan 7
23196 Malmö,
SWEDEN
info@Molderm.com

0373

Italfarmacia s.r.l.
Via di Tor Sapienza
7 00155 - Roma, ITALY

Symbols

-  Do not use if package is damaged
-  Warnings
-  Store between 2°C and 30°C
-  Do no reuse
-  Store away from light
-  Do not resterilize
-  Read the instructions for use
-  Sterlized using steam heat
-  Do not disperse into environment
-  Keep in a dry place

The SSCP (summary relating to safety and clinical performance) is available on the website of the European database on medical devices Eudamed (<https://ec.europa.eu/tools/eudamed>), where it is linked to the UDI-DI code of base 805571533ITF00578abdYB

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